CAMP[®]
SCANDINAVIA

AVLASTNING

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 **OPTIMA**[®]
molliter

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Franco Salvatelli R&D

Susanna Salvatelli C.E.O.

Alberto Salvatelli Product Manager

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Emerging Treatments and Technologies
ORIGINAL ARTICLE

An Off-the-Shelf Instant Contact Casting Device for the Management of Diabetic Foot Ulcers

A randomized prospective trial versus traditional fiberglass cast

ALBERTO PIAGGIO, MD
SEVIL MANGELDING, DPM
LORETTA RIZZI, MD
FRANCISCA PALMIERI, MD
ANNE TESCHKE, MP

LIAURO AMBRONI NORI, MD
ELISA LEPIATRI, DPM
VINCENZO SOBRI, DPM
IRASIA THOMAS, DPM
STEANO DEL PLATO, MEL, PDM

OBJECTIVE: This study was designed to test the safety, effectiveness, and costs of off-loading with a novel, off-the-shell removable device in the management of diabetic foot ulceration (DFUs).

RESEARCH DESIGN AND METHODS: We prospectively evaluated off-loading of neuropathic plantar ulcers in 40 diabetic inpatients attending our diabetic foot clinic and compared healing rates at the 12-week follow-up, number and severity of adverse events, healing time, costs and acceptability of the device, and patient satisfaction between those randomly assigned to instant contact casting (TCC; group A) or to the Optima-Dub walker (group B). Depth of infected ulcers were evaluated.

RESULTS: No difference between groups A and B was observed in healing rates at 12 weeks (35 vs. 37%), healing time (6.5 ± 4.4 vs. 8.7 ± 3.4 weeks), and number of adverse events (in seven patients). Treatment was significantly less expensive in group B, which showed a mean reduction of cost of 78% compared with group A ($P < 0.001$). Fractures were more favorable in group B, with a reduction of 77 and 58% of the time required for application and removal of the devices, respectively ($P < 0.001$). Patients' satisfaction with the treatment was higher in group B ($P < 0.01$).

CONCLUSIONS: The Optima-Dub walker is as safe and effective as TCC in the management of DFUs, but its lower costs and better applicability may be of help in spreading the practice of off-loading among the centers that manage the diabetic foot.

Diabetes Care 30:389–390, 2007

Foot & Ankle International

Comparison of Removable and Irremovable Walking Boot to Total Contact Casting in Offloading the Neuropathic Diabetic Foot Ulceration

Alberto Piaggi, MD¹, Chiara Goretti, MD¹, Elisabetta Iacopi, MD¹, Giacomo Clerici, MD², Fabio Romagnoli, MD³, Fabrizia Toscanella, MD⁴, and Cristiana Vermigli, MD⁵

Abstract

Background: Despite its efficacy in healing neuropathic diabetic foot ulcers (DFUs), total contact cast (TCC) is often underused because of technical limitations and poor patient acceptance. We compared TCC to irremovable and removable commercially available walking boots for DFU offloading.

Methods: We prospectively studied 60 patients with DFUs, randomly assigned to 3 different offloading modalities: TCC (group A), walking boot rendered irremovable (i-RWD; group B), and removable walking boot (RWD; group C). Patients were followed up weekly for 90 days or up to complete re-epithelialization; ulcer survival, healing time, and ulcer size reduction (USR) were considered for efficacy, whereas number of adverse events was considered for safety. Patients' acceptance and costs were also evaluated.

Results: Mean healing time in the 3 groups did not differ ($P = .5579$), and survival analysis showed no difference between the groups (logrank test $P = .8270$). USR from baseline to the end of follow-up was significant ($P < .01$) in all groups without differences between the groups. Seven patients in group A (35%), 2 in group B (10%), and 1 in group C (5%) (Fisher exact test $P = .0436$ group A vs group C) reported nonsevere adverse events. Patients' acceptance and costs were significantly better in group C ($P < .05$).

Conclusions: Our results suggest that a walking boot was as effective and safe as TCC in offloading the neuropathic DFUs, irrespective of removability. The better acceptability and lesser costs of a removable device may actually extend the possibilities of providing adequate offloading.

Level of Evidence: Level II, prospective comparative study.

CLINICAL CASES AND REFERENCES:

The International Journal of Lower Extremity Wounds
The Use of a Dermal Substitute (Integra) to Preserve Maximal Foot Length in a Diabetic Foot Wound
INT J LOW EXTREMITY WOUNDS 2006; 12: 200–201. © 2006 Blackwell Publishing Ltd

J BJS
ORTHOPEDIC SURGERY AND RESEARCH
The Journal of Bone & Joint Surgery – American Volume

WoundPedia™
Evidence Informed Practice

Wounds
INTERNATIONAL

THE AMERICAN
JOURNAL OF
MEDICINE®

AMERICAN ORTHOPAEDIC
FOOT & ANKLE SOCIETY
DISSEMINATION • PROFESSIONAL • PRACTICE

Biomechanics
For Motion in Daily Activities and Medicine

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DiabetesCare

Diabetes Care N°30 03/2007, p. 586-590

An Off-the-Shelf Instant Total Contact Casting Device for the management of Diabetic Foot Ulceration 15/2016

FAI FOOT & ANCLE INTERNATIONAL

Foot & Ankle International Vol 37, Issue 8, 2016

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