
elife


Application Instruction

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

Version EN (English)

MWKFS06 Solid Back Walker

- ① Place foot and shin inside foam liner. Then, fasten liner snugly.

Note:

You may put pad kits between walker frame and liner or inside liner to fill inner space on demand.

- ② Place wrapped foot and shin into walker frame. Then, place the heel against the posterior panel of the walker frame.
- ③ Place front panel over shin and ankle.
- ④ Fasten and tighten ankle, foot and shin strap in order snugly and comfortably.

Size Information

Length of inner sole.

Size	CM
S	25
M	28
L	30
XL	32

Washing instruction

- Remove straps, pads or liners and close all hook and loop closures before washing.
- Hand wash with mild detergent and cold water in gentle cycle.
- Foam liner is able to be machine washed at temperature not exceeding 40°C.
- Do not dry clean or tumble clean and do not use any bleach.
- Dry in the air current and not expose this brace to any heat source.

Caution

If allergic skin rashes or discomfort increases, please consult your physicians.
Do not use more than 24 hours of continuous.

Indication

- Grade II / III ankle sprain
- Stable fracture of the distal fibula, distal tibia, talus, calcaneus and metatarsal bones
- Post foot injury

Contraindication

Paraesthesia
Severe varicose vein
Skin disorders or injuries in the relevant area
Impaired lymphatic drainage

Recommendation

These products are to be used only by those which have been medical specialist. In order to dispose of the product, comply with the legal norms of your country. If any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Classification

Class I medical device, non sterile, no measurement function. This product fulfils the Medical Device Regulation (EU) 2017/745.

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e-life International Co., Ltd.
MDSS GmbH