Ortesis de cadera  
جهاز تقويم الورك  
Kýčelní ortéza  
Hüft-orthese  
Hofteortose  
Κηδεμονας του ισχιου  
Hip orthosis  
Lonkkaortoosi  
Orthèse de hanche  
Ortoza za kuk  
Csípőrögzítő ortézis  
Ortesi per anca  
Heuporthese  
Hofteortose  
Orteza stawu biodrowego  
Ortôtese da anca  
Ortežã de șold  
Ортез для тазобедренного сустава  
Ortéza bedrového kĺbu  
Höftortos

PSL699

Instrucciones de adaptación · تعليمات المستخدم · Pokyny K Nastavení · Anleitung ·  
Tilpasning instruktion · Οδηγίες εφαρμογής · Fitting instructions · Sovitusohjeet ·  
Instructions d'adaptation · Upute za postavljanje · Rögzítési utasítások · Istruzioni  
per l'adattamento · Montage-instructies · Tilpassingsinstruksjon · Instrukcja obsługi  
· Instruções de adaptação · Instructiuni de reglare · Инструкции по адаптации ·  
Návod na použitie · Utprovningsinstruktion

# ES

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## ENG

### Description and intended use:

Modular, one-size-fits-all hip stabilising orthosis that can be used for either the left or right hip. Composed of a pelvic band and thigh support connected by a joint that can be adjusted to provide from 0° to 90° of flexion-extension in intervals of 15°, and 0°, 15° or 30° of abduction.

### Indications:

Post hip surgery. Replacement or repair of total hip prosthesis. Hip dislocation. Fractures associated with hip replacement. Instability of the hip joint. Osteoporosis. Loosening of hip endoprostheses. Moderate to severe coxarthrosis.

### Warnings:

- Read these instructions carefully before using the orthosis. If you have any questions, consult your doctor or the shop that sold you the product.
- For the orthosis to serve its purpose, it must be correctly adjusted to the patient's morphology. The first fitting and adjustment must be carried out by qualified personnel (doctor, orthopaedic technician, etc.)
- Always follow these general instructions of use and any particular guidelines indicated by your general practitioner.
- Only a doctor is qualified to prescribe and determine the duration of treatment and carry out subsequent monitoring.
- If you notice any side effects, skin conditions or sensitivity, consult your doctor immediately.
- Your doctor might consider the use of topical creams in conjunction with the orthosis.
- In areas where pressure is applied, the skin should not be damaged or hypersensitive.
- Although the orthosis is not for single use, it should only be used by one patient. The product should not be reused by other patients.
- For the orthosis to function correctly, its components must be regularly checked to ensure they are in good condition. The medical staff responsible for monitoring the treatment may recommend the product to the patient or suggest that it be replaced if it is damaged or worn.
- To ensure they function correctly, clean the hook-and-loop fastenings regularly, removing any material that may have stuck to them.
- Some of the materials that make up the product are not flame retardant. Therefore, its use is not recommended in environments with excessive heat, fire or radiation.
- Clean the product regularly and maintain proper hygiene.
- Ask your doctor if you can remove the orthosis when bathing or showering. If you cannot, make sure you cover it or keep it away from the water.
- Do not let children play with this device.

### Fitting and adjustment instructions:

1. Adjust the pelvic band to fit the patient as shown in fig. 9 before fitting the orthosis.
2. Position and adjust the height of the joint.

2.1 Firstly, attach the upper shell of the joint to the pelvic band on the side of the hip to be treated using the hook-and-loop self-fasteners. Important: ensure the hook-and-loop self-fasteners adhere along the entire width of the band. Note: The orthosis comes preassembled for the right side.

Once the upper shell is in place, adjust the height of the joint by aligning the adjustment wheel with the height of the trochanter. To do this, use the Allen key to loosen the screws that attach the joint to the shell and move it until the desired position is achieved. After adjusting, retighten the screws.

If necessary, to improve the fit to the patient's anatomy, the inclination angle of the joint's upper bar can be adjusted. To do this, use an Allen key to loosen the upper side screw (A) (see fig. 1) without removing it completely, then adjust the inclination

to achieve an optimal fit and retighten the screw to secure. You can choose between 5 different positions.

2.2 Secondly, adjust the length in the lower shell area. To do this, loosen the screws that attach the joint to the lower shell and move it until the condylar support is aligned with the condyle in the patient's inner knee. Once the height has been adjusted, retighten the screws and fix the straps, as shown in fig. 2. Important: depending on the hip to be treated, the straps must be altered so that they are always in the position shown in fig. 2, i.e., the large buckles must always be on the front of the thigh. See fig. 3.

If necessary, more height adjustment can be achieved by changing the position of the lower shell's screws, as shown in fig. 4.

If internal rotation over the patient's knee is required, pull the straps behind the knee as shown in fig. 5. If necessary, trim the excess strap lengths, as shown in fig. 6. To do this, remove the crocodile fasteners, trim the excess strap and reattach

Once the entire orthosis has been correctly adjusted, pass the pelvic band tensioning straps between the shell and the joint as shown in fig. 7, tighten as required and use the self-fastening hook-and-loop strips to secure them to the belt. If the straps are too long, remove the crocodile fasteners, trim the excess, then reattach and adjust.

#### **Abduction adjustment:**

To adjust the degree of abduction, loosen the lower side screw (B) (fig. 1) without removing it completely. Next, select the desired angle ( $0^{\circ}$ ,  $15^{\circ}$ ,  $30^{\circ}$ ) and tighten the screw. For convenience, this adjustment should be made before fitting the orthosis to the patient.

#### **Flexion-extension adjustment:**

To adjust the flexion-extension, use an Allen key to loosen the central screws (C) (see fig. 1) until you are able to rotate the cogwheels (D) (see fig. 1) without releasing them completely. To make it easier to rotate the cogwheels, after loosening the screws, insert the tip of the Allen key into the hole marked by the arrows on the cogwheels, as shown in fig. 8. Then adjust the degrees of flexion and extension according to the patient's requirements. For convenience, this adjustment should be made before fitting the orthosis to the patient. After the adjustment has been made, re-tighten the screws.

#### **Contraindications:**

Not indicated for use on open wounds with swelling, reddening, or heat buildup.

#### **Product maintenance - washing instructions:**

To clean, close the hook-and-loop self-fasteners and wash with a neutral detergent (max.  $30^{\circ}\text{C}$ ). Do not use bleach. Do not expose to direct heat sources such as stoves, heaters, radiators, direct sunlight, etc.

To clean the metal structure and plastic shells, use a damp cloth and then wipe dry.

#### **Note:**

**Any serious incident involving the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.**

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PRIM, S.A.