








Aspen[®]

MEDICAL PRODUCTS

ASPEN[®] CERVICAL COLLAR

CERVICAL COLLAR
ZERVIKALORTHESE
CUELLO CERVICAL
COLLIER CERVICAL
NEKBRACE
KRCŃNÍ LÍMEC
HALSHVIRVEL KRAVE
KRCŃNÍ LÍMEC
CERVIKAL KRAGE
SZYJNY KOLNIERZ ORTOPEDYCZNY
KAULATUKI
NACKKRAGE
ШЕЙНЫЙ КОРСЕТ
COLAR CERVICAL
頸托

ENGLISH
1. **SIZING**
Maintain head in desired position. Place the sizing guide against head at the highest point of the shoulder muscle. Draw an imaginary line from chin bottom to the sizing guide then select collar size.
2. **PRE-FORM**
Roll up the back-panel.
3. **POSITION FRONT**
A. Flare sides of the front panel outward. Place chin piece directly under the chin. Generally, the chin should not extend beyond the edge of the plastic.
B. Push the sides of the front panel up over shoulder muscles around the neck.
4. **ATTACH FRONT**
A. While holding the front panel with one hand, center the back panel and attach both sides to the front.
B. While holding the front panel in place with one hand, pull outward then secure the back panel straps to the front of the collar. Repeat as needed for a proper fit.
5. **ADJUSTMENTS**
Tighten the Occipital Support Strap™ by squeezing the back panel to the neck and pulling until all slack is removed.
NOTE: It is not necessary to tighten the Occipital Support Strap™ on supine patients.
6. **PROPER FIT**
A. The patient's chin should be flush with the end of the collar chin piece. The inner trach bar should not be touching the airway. If it does, refit with the next taller size.
B. All slack should be removed from the collar back. The back panel should be centered. From the front, the back straps should be symmetrical.

INTENDED USE
To provide motion restriction of the cervical spine.
INDICATIONS
Decompressive procedures, Fracture management, Postsurgical stabilization, Whiplash, Sprains/strains, Cervicogenic headache, Cervical disc syndrome, Musculoskeletal pain, Radiculopathy.
CONTRAINDICATIONS
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.
PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
Hand wash only at 30° C. Wash with mild detergent; Air dry; Do not use bleach or other chemicals.
MATERIALS
HDPE, Nylon, Polyurethane, Cotton.
NOT MADE WITH NATURAL RUBBER LATEX.
DISPOSAL
Product may be safely disposed in accordance with local laws.
REPORTING
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).
LIMITED WARRANTY
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
(a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
(b) The product is altered;
(c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
(d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.
THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OR DEFECT IN THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES OF THE PRODUCT, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.
This warranty gives you specific legal rights and you may have other rights which vary from location to location.
A listing of authorized distributors in your area is available upon request.

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