Orthelligent[®] by OPEd VISION

BECAUSE EVERY STEP COUNTS.

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Istruzioni per l'uso

OPEd
Keeps you going.

DEUTSCH

ENGLISH

FRANÇAIS

ITALIANO

DEUTSCH

Willkommen

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Orthopaedic Appliances Pty Ltd
OAPL an OPED company
26-32 Clayton Road
Clayton, Vic 3168
Australia
Phone +61 1300 866 275
Fax +61 1300 627 539
info@oapl.com.au
www.oapl.com.au

OPED Medical, Inc.
444 Tom White Road
Braselton, GA 30517
USA
Phone +1 (770) 945 – 0150
Fax +1 (770) 945 – 0370
info@opedmedical.com
www.opedmedical.co

6 Symbolbeschreibung

	EN	DE	FR	IT
	Medical Device	Medizinprodukt	Dispositif médical	Dispositivo medico
	Manufacturer	Hersteller	Fabricant	Produttore
	Unique Device Identification	Eindeutige Geräteerkennung	Identifiant de dispositif unique	Identificatore univoco del dispositivo
	Read the operating instructions	Gebrauchsanweisung beachten	Consulter le mode d'emploi	Osservare le istruzioni per l'uso
	Swiss authorised representative	Schweizer Bevollmächtigter	Mandataire suisse	Rappresentante autorizzato svizzero

ENGLISH

Welcome to

Orthelligent® VISION is a tablet application (app) for analyzing movements such as walking, including a written evaluation. The app supports the measurement and documentation of movements, calculates recognized parameters and provides feedback based on the current state of science and medical guidelines.

Please read these instructions carefully before starting Orthelligent® VISION. The current version of the app supports you in carrying out and documenting movement analyses using the tablet camera. The explanations in the following chapters only cover the functions and areas of application of the current version of the app.

We have summarized important information for a smooth start here. If any questions are not answered here, please contact us by email at support.dhs@oped.de.

An introductory video is available on our website www.o-dhs.com or on [Orthelligent®- YouTube](#).

1 Information on Orthelligent® VISION

1.1 Intended use

Software for the documentation of dynamic and static motion sequences for professional users.

1.2 Usage environment

Orthelligent® VISION is designed for use in the healthcare sector, taking into account the safety instructions described in the app and the instructions for use. The measurement can be carried out both indoors and outdoors. An area of approx. 5 x 5 m should be available.

1.3 User and profile

This medical device has been developed for use by healthcare professionals such as doctors, physiotherapists, technicians or similar professional groups.

No special training is required, the app is self-explanatory.

1.4 Technical requirements for using the Orthelligent® VISION app

The app can only be installed on a tablet with an iOS or Android operating system.

Minimum requirements for the tablet:

- Operating systems: Apple iPad: iOS version 11 or higher;
Android version 7.0 or higher
- Internet connection (after installing the app, operation in offline mode is possible to a limited extent)
- At least 1.5 GB RAM working memory
- At least 100 MB usable free disk space
- Please protect the app from access by unauthorized persons.

1.5 Warnings

- A movement analysis by Orthelligent® VISION cannot diagnose a disease and is not a substitute for medical advice.
- Orthelligent® VISION is not an alarm system for emergencies. If individual measured values appear unusually high or low, they should definitely be clarified by a doctor.
- The questions about the patients must be answered carefully and correctly. This serves to ensure the safety of the measurements. It is essential that the patient's height is stated. Incorrect information can lead to incorrect measurements and analyses.

1.6 Indications

Orthelligent®VISION can be used to analyze and assess movements. Measurement is possible both with aids (e.g. with an orthosis) and without.

1.7 Contraindications

There are no contraindications. However, the general state of health should permit a movement analysis. Furthermore, the patient should be able to walk about 5 m alone.

Side effects: There are no known side effects.

2 Registering and starting the app

2.1 Registration

An e-mail address that can be accessed at any time is required for registration. Please make sure that no unauthorized persons have access to it. To create your account, visit the O-DHS webshop and register with your e-mail address.

Register your user account at:

<https://webshop.o-dhs.com>

2.2 Installation in the app

Install the app via Google Play (Android) or the AppStore (Apple). Start the app to log in. To log in, you will need your e-mail address (username) and password, which were assigned when you created your account. Then follow the further instructions in the app. If you do not yet have a user account, you can register at any time via the webshop and create an account.

2.3 Registration (login) in the app

After registering in the app with your e-mail address and password (see 2.1 Registering in the app), you can log in to the app at any time using your access data.

2.4 Logout

In addition to the account button with the user's initials, which can be found at the top right, you can log out manually at any time. To protect against unauthorized access, the logout occurs automatically after 10 minutes without input.


Fig. 1: App dashboard

2.5 Deactivation of the account

To deactivate your account, send an email to support@o-dhs.de.

2.6 Uninstalling the app

Android tablet: Go to your home screen or open the app drawer where all your apps are displayed. Find the app you want to uninstall. Tap and hold the app icon until a menu or option appears. To uninstall, either drag the app icon to the top of the screen, where "Uninstall" or a trash can icon appears, or select the "Uninstall" option in the pop-up menu. Then confirm the uninstallation.

iOS tablet (iPad): Go to your home screen, where all your apps are displayed, and find the app you want to uninstall. Tap and hold the app icon until all the icons start to wiggle and a small "X" appears in the top corner of the icon. Tap on the "X" of the Orthelligent®VISION app and confirm the uninstallation. End the "wiggle mode" by pressing the home button or tapping "Done" (on newer iPads without a home button).

3 Home page and dashboard

The page shown in Fig. 2 appears at the beginning.


Fig. 2: Positioning a patient

- To create a new patient, click on the "Add patient" button.
- The patient list can be viewed via the "Patient" icon on the left-hand side.
- Profile: information about the user account can be viewed here
- Help: problems or errors can be reported here. Contact details are provided.
- About us: Here you can find the privacy policy, the imprint and further details.
- The dashboard can be accessed at any time via the **VISION** icon at the top left.

4 Carrying out movement analyses

Movements are multifaceted, and even small discrepancies often have a major impact on well-being, regeneration and injury prevention. The markerless movement analysis of Orthelligent® VISION enables a precise analysis of movements including a standardized evaluation.

Flow chart:

- > **Create or call up patient profile**
- > **Create analysis**
- > **Start analysis**
- > **Record video**
- > **Check and confirm video**
- > **Upload video for analysis**
- > **Receive video with markers**
- > **Check and confirm data**
- > **Receive and share report**

4.1 Creation of the patient profile

To be able to assign the patient after the measurement, the patient's profile must be created. Press the ⊕ icon in the patient list or the "Add patient" button to create another profile (see 4.1.)

The following information is required to create a profile:

- Patient number (this can be freely selected, according to the existing number range) *
- Surname, first name and gender
- Date of birth *
- Size, must be specified (exact specification is relevant for further calculations)*
- Weight (can be deselected, see Fig. 3)
- Type of physical activity
- Information on types of sport practiced and frequency of sporting activity
- Injury incl. date and type

All fields marked with * must be completed for a reliable analysis.

After entering the data and checking the information in the last view, click the "Add patient" button to complete the process.

4.2 Measurement

General requirements

- The person's head and feet must be in the picture throughout. Only one person may be seen in the video.
- Minimum body height of 90 cm required.
- The patient should be able to walk 10 steps straight ahead alone.
- Patient should wear tight-fitting clothing without patterns (no dots/stripes). Skirts and very dark clothing can falsify results.
- Align the camera to the center of the body.
- Hold the camera at least 2 m away, stabilize the camera.
- Adequate lighting without strong background lighting.
- The silhouette on the tablet serves as an orientation for the optimum position of the test subject in the video. If the silhouette lights up green, the tablet is in the optimum position.
- Alternatively, a treadmill/walking pad can also be used for the measurement. It is important that the head and feet are clearly visible.
- Always have the analysis carried out by a qualified specialist. Do not place the tablet on the floor to avoid distortion.
- Up to 10 measurements can be recorded in offline mode. The results are made available as soon as the tablet is back online.

Specific requirements for frontal analysis and staircase analysis:

- Ideally in a room measuring 2 m x 5 m.
- At least 5 steps per analysis in one run. Ideally, walk towards the camera twice, 5 steps each time. Walking away from the camera is not analyzed.
- Switch off the rotation lock in the iPad. For older iPads, check the side switch and select the side switch function in "Settings" > "General".

Specific requirements for lateral analysis:

- Ideally in a room measuring 5 m x 5 m.
- At least 6 steps per analysis in one run.

After a successful measurement, follow the instructions in the app to upload and calculate the video using the upload icon.

Measurement procedure

Open the patient profile and start a new measurement. At the beginning, select which measurement is to be performed. (Fig. 3).

The following options are currently available:

Lateral, frontal, stairs


Fig. 3

Now click on the \oplus icon at the top right to create a new flow for your measurement. You can give it any name you like.

After clicking on the "Save and exit" button, the analysis can be started. To do this, simply open the "Open flow" and start the "Start scan".

Existing videos can also be imported into our app and then analyzed if they meet the above requirements. To do this, simply click on the \oplus icon after starting the analysis and select the file from the device. (Fig. 4)


Fig. 4

However, it is important here that the video meets the general requirements of our app, which are listed below. These instructions are also relevant for the regular measurement within the app:

Completion of the measurement:

By clicking on the upload icon, the video is then analyzed and the virtual markers are placed on the patient. (Fig. 5)

Fig. 5

4.3 Report and evaluation

The automated report is at the heart of markerless motion analysis with Orthelligent®VISION.

4.3.1 Creation and confirmation

In order to receive the report of the markerless motion analysis with Orthelligent®VISION, the joint points set by the AI must first be checked to ensure that the person in the video has been correctly recognized. To do this, simply click on "Check results" after the upload process as soon as an orange symbol appears to display the video with the virtual markers. If a red symbol is displayed, the measurement has unfortunately not worked and you will have to carry it out again, pay attention to the error message for the video.

The report can be downloaded and shared using the three dots at the top right (Fig. 6).


Fig. 6

The report can be downloaded to the device and saved or sent (Fig. 7)


Fig. 7

4.3.2 Structure and content

The report is divided into four parts. At the beginning, general parameters are introduced, followed by an overview of all gait parameters. This is followed by a detailed evaluation of the 8 gait phases and a presentation of the angle curves.

4.3.3 Parameters

The following parameters are listed in the report:

Number of steps: The number of steps the patient took during the analysis.

Gait symmetry: This is calculated in the lateral measurement via the symmetry of the left and right stride lengths. The higher the value, the better.

In the frontal analysis, this value represents the average stride time comparing the left and right leg.

Gait speed: shows how fast the patient walked during the analysis.

Double support: Indicates how often the patient stood with both feet on the floor during the analysis.

Gait variability left/right: Indicates the fluctuation of the stride length of the left/right leg during the analysis.

Cadence: indicates the total number of steps taken within one minute during the analysis.

Step length left/right: indicates the average distance between the point of first contact of the left foot and the point of first contact of the right foot during the analysis.

Step time left/right: shows the average time that elapses during the analysis from the first contact of one foot to the first contact of the other foot.

Stride width left/right: The stride width describes the lateral distance between the heels of both feet during walking. It measures the width of the stride pattern and can indicate gait stability and balance.

4.4 Evaluation further information

In addition to the parameters, the 8 gait phases are examined in detail. The corresponding angles in the hip, knee and ankle joints are displayed for each phase and shown in relation to literature values.

The first page of the report summarizes the measured values in a short text and also provides space for notes on the measurement.


Fig. 8

The second page shows the space-time parameters. The values are classified in the context of the reference values.


Fig. 9

The third page shows the eight gait phases according to Götz-Neumann, and the specific measured values for each gait phase are illustrated for the hip, knee and ankle in both the left and right leg. The marked area within the bars represents the reference value.


Fig. 10

On pages four and five, the measured values are again assigned to the individual gait phases and explained in a summarizing text.


Fig. 11


Fig. 12

The last page of the report shows the flexion and extension of the respective joints and displays the values of the left and right leg in comparison to the reference range in a diagram.


Fig. 13

The report can be downloaded from the app using the "Export" button.


Fig. 14

4.5 Comparison of two measurements

As soon as you have two reports for one type of analysis (lateral, frontal, stairs), you can compare the measurements. To do this, click on the desired analysis, then on the three dots in the menu field and then select "Compare".


Fig. 15

You can freely select the reference measurement and change it at any time. Then give the comparison a suitable name, such as "Comparison 1" or a name of your choice.


Fig. 16

5 Contact and further information

We will be happy to answer any further questions at dhs.support@oped.de

OPED GmbH
Medical Park 1
83626 Valley/Oberlindern, Germany

Phone +49 8024 60818-210

Fax +49 8024 60818-299

mail@oped.de

Data protection


Information about the processing of your personal data can be found at:

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According to MDR (EU) 2017/745 Annex VIII, Orthelligent®VISION is a Class I medical device.

In the event of complaints in connection with the product, please contact the manufacturer directly. Only serious incidents that may lead to a significant deterioration in health or death must be reported to the manufacturer and the competent authority of the member state. Serious incidents are defined in Article 2 No. 65 of Regulation (EU) 2017/745 (MDR).

International contacts:

OPED AG
Hinterbergstrasse 26
6312 Steinhausen, Switzerland
Phone +41 41 748 53 35

contact@oped.ch www.oped.ch

Orthopaedic Appliances Pty Ltd
OAPL an OPED company
26-32 Clayton Road
Clayton, Vic 3168
Australia
Phone +61 1300 866 275
Fax +61 1300 627 539
info@oapl.com.au
www.oapl.com.au

OPED Medical, Inc.
444 Tom White Road
Braselton, GA 30517
USA
Phone +1 (770) 945 - 0150
Fax +1 (770) 945 - 0370
info@opedmedical.com
www.opedmedical.co

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The screenshot shows a web-based report interface for a lateral gait analysis. At the top, there are navigation icons (back, forward, search) and an 'Exportieren' button. The main content area is titled 'Report: Ganganalyse (lateral)' and includes the ODHS logo (INNOVATION BY OPED). Below the title, patient and analysis details are listed in a grid format. A summary section titled 'Zusammenfassung' contains a paragraph of text describing the patient's gait characteristics. There are also sections for 'Notizen' (notes) and 'Bericht erstellt für:' (report created for: OPED GmbH) with a signature line. The footer contains a small ID '1b5e14f', a reference to a book by Götz-Neumann, and the page number 'Seite 1 von 6'.

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