THE BOSTON T

The Boston T is the perfect mate to wear under your Boston Brace. Made in the USA especially for Allard Support UK Ltd by Knit-Rite Inc. The Boston T is a round neck protective body sock that is superior to a regular t-shirt for cushioning. Stretches more than 150% to ensure a wrinkle-free fit under vour brace.

SIZE GUIDE:

Size	Width	Circumference	Length
Small	6 inches / 152mm	16-24 inches / 406-610mm	19 inches / 483mm
Medium	8 inches / 203mm	24-32 inches / 610-813mm	21 inches / 533mm
Large	10 inches / 254mm	32-40 inches / 813-102mm	28 inches / 711mm
X-Large	12 inches / 305mm	40-48 inches / 102-122mm	28 inches / 711mm

MATERIAL CONTENT:

Item Code Prefix	Description	Ply	Material Content
1BX1BV	lonic+™, Single Flap	1-ply	49% Polyester, 46.5% Nylon; 5% lonic+ [™] , 2% LYCRA® IONIC+ [™] fibres help transport moisture away from the skin and prevent odour in the garment. It also has thermal and anti-static properties to keep you cool and comfortable.
1BP2BV	Standard, Single Flap	2-ply	95.6% COOLMAX®, 4.4% LYCRA® COOLMAX® includes an effective fibre-based moisture management system. The system can move perspiration away from the body, and through the fabric, where it can evaporate quickly, allowing the wearer to feel cooler and more comfortable.
1BXBBD	Ionic+™, Double Flap	2-ply	95.6% Polyester / IONIC+™ blend (See material properties above), 4.4% LYCRA®

COOLMAX® and LYCRA® are registered trademarks of The LYCRA Company. Ionic+™ is a trademark of Noble Biomaterials, Inc.







WASHING INSTRUCTIONS:

Machine wash warm and tumble dry low. Only non-chlorine bleach (when needed).



Sarl. 11 rue Emile Zola. BP 2332 38033 Grenoble Cedex 2. France

#0053 2208



www.allarduk.co.uk