

VRTX Configurations

VRTX CTLSO Options

LSO + CTLSO Upgrade = CTLSO
 TLSO + CTLSO Upgrade = CTLSO

VRTX Lower Spine Options


LSO
 TLSO

Size	Circumference
X-SMALL	66 - 91.5 cm
STANDARD	78.7 - 152.4 cm
X-LARGE	129.5 - 177.8 cm

VRTX Upper Spine Options


48068 Collar, 48065 Collar or 48074 Collar
 CTO (Uses Multipost Collar 48065)
 CTO 4-post (Uses Multipost Collar 48065)

VRTX COLLAR




Part Description: VRTX Collar 48068
Type: Universal
Height Settings: 1 - 8 (21 total adjustments)
Height Range (Sternal notch to chin): 2,54 cm to 12,70 cm
Circumference: 33,02 - 53,34 cm


Anterior Panel Height VRTX LSO



Posterior Panel Height VRTX LSO and TLSO



Anterior Heights VRTX TLSO



Anterior Heights VRTX CTLSO

Posterior Heights VRTX CTLSO

